

RNA module for Countable PCR (KT0010)

For Research Use Only. Not for use in diagnostic procedures.

Description

RNA Module for Countable PCR Kit is designed to perform reverse transcription and PCR amplification from RNA templates within a single reaction. Optimized for use with the Countable 4X PCR Mix (KT0004) and the Countable PCR platform (IM0001), this module enables a streamlined one-step workflow that improves efficiency, specificity, and sensitivity for precise RNA quantification, while ensuring RNA integrity throughout the entire workflow using RNase Inhibitor. The kit is compatible with hydrolysis probe-based assays, delivering clean signal profiles and uniform amplification across a wide dynamic range.

Kit Contents

Component	Catalog #	Quantity
Reverse Transcriptase	PR0019	1
RNase Inhibitor	PR0018	1

Supports up to 192 Countable PCR reactions per kit.

Storage and Stability

- Store all components at -25°C to -15°C
- Avoid repeated freeze-thaw cycles
- Shelf life: 2 years from date of manufacture

Handling Instructions

- Thaw components on ice or ice-block and mix by gentle tapping. Centrifuge briefly to collect contents at the bottom of the tube
- Prepare the RNA sample before setting up the reaction mix and keep both of them on ice
- Suggested input quantities of total RNA are up to 1 µg per reaction

Workflow Overview

See Countable PCR Reaction Preparation User Guide

Safety Information

This kit does not contain a reportable quantity for hazardous materials. Use in accordance with standard lab safety practices.

Additional Information

Probe Additive Recommendation

For hydrolysis probe assays, Countable Labs strongly recommends the use of a Probe Additive (PA) to reduce background fluorescence and enhance signal-to-noise ratio. The Probe Additive is an oligonucleotide designed to improve quenching efficiency in single-molecule PCR reactions, where residual unquenched fluorophores can otherwise contribute to false signal detection.

Addition of Magnesium Chloride (MgCl₂)

MgCl₂ is an essential cofactor for enzymatic activity. The Countable 4X PCR Mix kit (KT0004) includes a separate 50 mM MgCl₂ solution. To enhance the signal in an RT-PCR reaction within the Countable system, add 1.5 µL of the 50 mM MgCl₂.

Recommended primer and probe concentrations are provided in the Countable PCR Reaction Preparation User Guide.

Legal Notice

All products from Countable Labs (including associated consumables and software) are for **Research Use Only**. They are not intended or approved for use in humans or animals, for diagnostic procedures, for clinical decision-making, or for any therapeutic or other medical purpose.

The products and/or their uses may be covered by one or more U.S. and/or foreign patents owned by or under license to Countable Labs. For patent information, visit: countablelabs.com/patents.

Countable Labs reserves all its intellectual property rights in and to its products, and associated consumables, software, and documentations. “Countable Labs”, “Countable PCR”, “Enumerix” are trademarks of Countable Labs, Inc. All other trademarks used herein are the property of their respective owners.

Unless otherwise provided in a separate written agreement executed by Countable Labs, Inc., use of a product from Countable Labs is subject to Countable Labs's applicable standard terms and conditions, available at: countablelabs.com/legal, which contain among other things, warranty disclaimers, liability limitations and use limitations, to the exclusion of any other terms. By opening or using the applicable product, you agree to be bound by those terms.

© 2026 Countable Labs, Inc. All Rights Reserved.