

Countable System Instructions for Use

For Research use only. Not for use in diagnostic procedures.

Countable Labs, Inc. 1810 Embarcadero Rd Ste 200 Palo Alto, CA 94303 USA countablelabs.com

Notices

This document is intended for the use of the Countable System. All products, including instruments, consumables, accessories, reagents, and services described herein are intended molecular biology applications and are FOR RESEARCH USE ONLY. They are NOT FOR USE IN DIAGNOSTIC PROCEDURES, or for the prevention and/or treatment of a disease.

Trademark and copyright

@2025, Countable Labs, Inc. All rights reserved.

Countable Labs, the Countable System, and the Countable Control Software are trademarks and/or registered trademarks of Countable Labs. All other trademarks are the property of their respective owners.

Duplication and/or reproduction of all or any portion of this document without the express written consent of Countable Labs is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. Countable Labs provides no warranty and hereby disclaims any and all warranties as to the use of any third party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product.

Countable Labs may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. The use of products described herein is subject to terms that have been agreed to in writing between Countable Labs, the user, and the user's respective organization.

Technical support

Your success is our success, so we are committed to supporting you at all times. Our Customer Success Team is made-up of scientists that are molecular biologists and Countable Labs product experts that are ready to answer any of your questions or help with any issue you may run into.

Contact us at success@countablelabs.com

Policy statement

We're always working to take our products to the next level so you'll have the best. Specification may be updated from time to time to reflect system improvements, but we'll keep you informed along the way. We also welcome your feedback on our guides so we can keep improving your experience.

Revision history

Document Revision	Revision Date	Description
IFUOO3 Rev 1.0	June 2025	Initial release

Table of contents

Introducing, the Countable System!	07
Welcome	07
Intended use of the Countable System	08
The Countable approach	08
Countable users	09
Preparing your lab	10
Overview	10
Space requirement	10
Electrical requirements	12
Environmental operating conditions	13
Safety and regulatory compliance	14
Proper use	14
Electrical safety	15
Mechanical safety	16
Instrument dimensions	16
Biological safety	17
Chemicals	18
Maintenance safety	18
Countable system compliance	19

General system description	20
System and workflow overview	20
The Countable instrument	20
The Countable computer	22
The Countable software	23
Countable consumables, accessories, and reagents	23

Operating the Countable System	25
Turning on your Countable System	25
Countable Control Software overview	27
Building a Sample Workbook	27
Defining Sample Workbook experimental parameters	30
Sample imaging	33
Viewing experiment reports	41
Countable Analysis Summary	43
Output data files	45
Refilling Countable Fluid	49
Emptying the waste bottle	52

Maintenance	53
Cleaning the Countable System	53
Cleaning carousel spills	54
Removing equipment from use for repair or disposal	54

Troubleshooting	55
Sending data to Countable Labs for troubleshooting or support	55
System initialization errors	55
Workbook is not accepted	56
Fluidics	57
Analysis summary issues and observations	58

CHAPTER 1

Introducing, the Countable System!

Welcome

Congratulations on bringing the Countable System to your lab! We're excited to welcome you to the next advancement of PCR technology, where enhanced precision and efficiency await you. To help you get the most out of your Countable System, we recommend taking a quick look at this Instructions For Use guide for system setup, operation, and maintenance information. For details on sample prep or assay-specific workflows, check out Countable PCR Reaction Preparation User Guide.

Conventions

Note	Note highlights handy tips to help you use your Countable System.
Warning	Warning calls your attention to situations where using the Countable System the wrong way could hurt you or someone else.
Caution	Caution alerts you to situations that might damage the Countable System or other equipment.

Intended use of the Countable System

The Countable System is intended for the absolute counting of target nucleic acids and genomic biomarkers. It is intended to be used by research and commercial laboratories to support their activities in research, development, and manufacturing. Applications included, but are not limited to, pathogen detection, targeted mutation detection, molecular residual disease testing, genotyping, copy number analysis, gene integration characterization, viral titration determination, methylation analysis, and treatment or therapy selection.

The Countable System is intended for Research Use Only (RUO).

The Countable approach

The Countable System is a high-resolution imaging and analysis platform that performs 3D light sheet fluorescence image of the Countable Matrix. It gives you precise spatial detection and counting of discrete fluorescent signals across the full reaction volume, making it the perfect solution for workflows that rely on single-molecule or compartmentalized fluorescence detection.

How does it work? Fluorescently-labeled molecules in your sample are immobilized in the Countable Matrix, a structured gel that forms in the Matrix Tube Strip. The Countable System images the tubes in the Matrix Tube Strip without any hands-on steps—giving you a no-transfer workflow from signal generation to image capture and analysis.

When your sample is imaged, the system collects a full-volume fluorescence map of each tube in the strip by scanning multiple focal plans. The Countable Control Software automatically generates signal detection, thresholding, and spatial compartment data, giving you high-resolution quantitative results.

Countable users

Task	Personnel	Required training and experience
Delivery	No special requirements	No special requirements
Installation and commissioning	Countable Labs service personnel	Trained Countable Labs personnel
Routine use (e.g. running protocols)	Laboratory technicians or equivalent	Personnel who are appropriately trained, or experienced and familiar with the use of computers and automation systems
Assay design and validation	Scientist or equivalent	Personnel who are appropriately trained or experienced and familiar with molecular biology techniques.
Cleaning and preventative maintenance	Laboratory technicians or equivalent	Personnel who are appropriately trained or experienced and familiar with the use of computers and automation systems
Service	Countable Labs service personnel	Trained and authorized Countable Labs personnel

CHAPTER 2

Preparing your lab

Overview

For a smooth and stress-free installation, review this chapter so your lab will be ready when a Countable Labs representative arrives to install your instrument. You should not try to install or move the instrument yourself.



Space requirement



Risk of personal injury and the Countable System damage

- For proper ventilation, maintain a minimum clearance of 4 in (10.2 cm) on the sides and rear of the Countable System.
- Don't cover any slits or openings that are essential for ventilation.

Instrument dimensions

Component	Length/Depth	Width	Height	Weight
Countable Instrument	24.8 in (63.0 cm)	19.7 in (50.1 cm)	13.6 in (34.5 cm) 22.8 in (57.9 cm) w/ door open	75.0 lbs (31.8 kg)
Desktop computer	16.5 in (42.0 cm)	6.8 in (17.3 cm)	14.5 in (36.9 cm)	18.7 lbs (8.5 kg)
Computer monitor	6.0 in (15.2 cm)	21.2 in (53.8 cm)	16.3 in (41.4 cm	7.3 lbs (3.3 kg)
Keyboard and mouse	7.4 in (18.8 cm)	17.6 in (44.7 cm)	0.8 in (2.0 cm)	1.1 lbs (0.5 kg)

Lab bench requirements

The recommended lab workbench is at least 30×72 (L" x W"), and the minimum table size is 30×60 (L" x W"). The lab workbench must be able to support 110 lbs (50 kg) to accommodate the weight of the Instrument and computer set up.

The imaging process in your Countable Instrument is sensitive to mechanical vibrations. It is very critical that your instrument is installed on a sturdy, wobble-free, vibration-free lab workbench for optimal image quality.

If you have control over the workbench being used, we recommend a Bench-Tek model # BT40CR-3072BS-PS, However, any existing lab workbench with similar footprint, quality, and sturdiness will also work just fine.

Instrument positioning requirements

You'll want to leave enough space around the side of your Countable Instrument for proper ventilation and so service personnel can access the instrument when needed. Minimum clearance requirements are shown below.

Side of instrument	Minimum clearance
Left side	At least 6.0 in (15.2 cm) clearance from the wall is required on the left side of the instrument, where the fan is mounted, for proper ventilation.
Right side and back side	4.0 in (10.2 cm) clearance is recommended for the right and back side of the instrument.
Top side	Shelves should be at least 27.0 in (68.9 cm) from bench surface.

Electrical requirements

Plug your Countable Instrument into a dedicated, grounded circuit capable of delivering 100–240 VAC, 50/60 Hz, 5 As.

Environmental operating conditions

Condition	Acceptable range
Installation site	For indoor laboratory use only.
Operating temperature	19—25 °C (66—77 °F).
Operating humidity	20—80% RH (non-condensing).
Altitude	<2000 m (6500 ft).
Air quality	Pollution Degree II or better.
Vibration	Don't place the instrument adjacent to strong vibration sources (e.g., centrifuges). Excessive vibration can negatively affect instrument performance.
Heat sources or direct sunlight	Don't place the instrument adjacent to heat sources or under direct sunlight. Temperature elevation inside the instrument due to heat sources or direct sunlight can negatively affect instrument performance.

CHAPTER 3

Safety and regulatory compliance

Before using your Countable System, review the following instructions and safety guidelines to avoid hurting yourself or damaging your instrument. The information in this chapter supplements, but does not replace, the safety requirements of your laboratory and/or institute.

Proper use	
	Risk of personal injury and the Countable System damage
Warning/	You should only operate the system after you've been trained.
Caution	 Only Countable Labs personnel should perform instrument service and maintenance.
<u>[1]</u>	Do NOT attempt to lift the Countable System yourself. The instrument should only be moved by Countable Labs personnel.
	Do NOT attempt to move the Countable System during sample imaging and analysis.

Risk of the Countab	le S	ystem	damage
---------------------	------	-------	--------

Caution

- Take care not to spill and liquid of any kind onto the Countable System. Damage caused by spills will void your warranty. In an emergency, turn off the instrument and unplug the power cord from the outlet.
- Only use Countable Labs validated consumables, accessories, and reagents with your Countable System. Damage caused by using other consumables may void your warranty.
- Don't put objects other than designated Matrix Tube Strips into the Countable System carousel.
- Don't stack or place objects on top of the Countable Instrument.

Electrical safety

General guidelines

- Connect the power cord to an outlet with proper grounding.
- Don't adjust and/or replace internal parts of the instrument yourself.
- **NEVER** operate the instrument with any covers or parts removed.
- If the instrument becomes electrically unsafe, stop others from using it and contact Countable Labs Support at success@countablelabs.com.

Don't operate the Countable System if:

- The instrument or power cord is or appears to be damaged.
- It has been stored in unfavorable conditions for a long time.
- Liquids have come into contact with electrical components on the instrument.

Warning Before powering on the Countable System for the first time, make sure you're using the correct supply voltage to avoid damaging the instrument electronics. Always connect the equipment to a grounded power outlet with the appropriate voltage and current rating.

Risk of the Countable System damage: Electrical hazard

Warning



- There are lethal voltages inside the instrument. When the instrument is connected to power, internal terminals could be live.
- Don't open covers and/or remove internal instrument parts when the instrument is connected to power, as this could expose you to live components.

Warning/ Caution

Risk of personal injury and the Countable System damage: Electrocution

- Don't open any panels on the Countable Instrument.
- Vou should only perform the maintenance steps described in this guide. Any other maintenance and/or repairs should only be carried out by Countable Labs personnel.

Mechanical safety

Warning/ Caution



Risk of personal injury and the Countable System damage

- Don't open any panels on the Countable Instrument when you're operating it.
- You should only perform the maintenance steps described in this guide. Any other maintenance and/or repairs should only be carried out by Countable Labs personnel.

Instrument dimensions

Safety standards

■ The Countable Instrument has been tested to IEC 60825-1:2024 & EN 60825-1:204/A11:2021.





Class 1 Laser Product

Risk of personal injury

- Caution Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure.
- You should only operate the system after you've been trained.
- The Countable Instrument is a CLASS 1 LASER PRODUCT. Never operate your Countable System with any cover panels removed and/or the system cover open.
- If the laser appears damaged or unsafe, turn off the instrument immediately and contact Countable Labs support.

Biological safety

Any specimens and reagents that you obtained from humans should be treated as potentially infectious. You should follow the standard laboratory practices for your institute and wear appropriate personal protective equipment when handling biological samples, as outlined in the CDC Biosafety in Microbiological and Biomedical Laboratories Manual (www.cdc.gov/labs/bmbl/index.html).

Samples

	Risk of personal injury: Samples containing infectious agents
Warning	Samples used with the Countable System may contain infectious agents. Handle all samples with caution and follow the required safety regulations.
	Always wear safety glasses, gloves, and a lab coat.
	You should only operate the system after you've been trained.
	Refer to relevant Material Safety Data Sheets to avoid being exposed to hazardous levels of infectious agents.

Personal Protective Equipment (PPE)



Chemicals

Warning Some Countable Consumables used with the Countable System may be hazardous. Handle reagents and chemicals according to their Safety Data Sheets. SDS' can be found on the Countable Labs website. Always wear safety glasses, gloves, and a lab coat.

Maintenance safety

Warning <u>()</u>	 Risk of personal injury and the Countable System damage The only maintenance you'll need to perform is described in this guide → See <u>"Chapter 6: Maintenance"</u> on page 58 for more information. Additional maintenance should only be performed by a Countable Labs personnel.
Warning	Risk of the Countable System damage
\bigwedge	Don't use bleach, solvents, reagents with acids, alkalies, abrasives, alcohol, or disinfectants to clean the Countable System unless otherwise instructed.

Countable System compliance

The Countable System is certified to the following standards:

EN 61010-1:2010/A1:2019 (Edition 3.1), IEC 61010-1:2010/AMD1:2016/COR1:2019, IEC 61010-2-081:2019, CSA C22.2 No. 61010-2-081:19

CHAPTER 4

General system description

System and workflow overview

Your Countable System includes the Countable Instrument, the Countable Computer, a monitor, a mouse and a keyboard. The Countable Computer comes ready to go with Windows 11 and the Countable Control Software already installed. You can tap into Countable technology by using your Countable System with Countable consumables, accessories, and reagents. Check out the Countable PCR Reaction Preparation User Guide for more information on using Countable consumables and reagents.

The Countable Instrument

Your Countable Instrument is a walk-away analytical instrument that takes care of the post-amplification steps for you—automatically detecting and measuring target molecules in your molecular biology assay so you can focus on the bigger picture. Matrix Tube Strips in the Countable System are imaged with lightsheet microscopy to give you high-resolution, 3D imaging in up to four channels of your choosing. This leads to rapid detection of single-molecule signaling across the entire 50 µL reaction.

Caution

Risk of error

Avoid bumping or causing vibration to your Countable System or the bench it is placed on, as it may interfere with system performance.



Labels on the Countable Instrument

Label description	Label image
Countable Fluid (lower): located on the work surface in front of Countable Fluid bottle position.	COUNTABLE FLUID
Countable Fluid (upper) : located on the vertical surface behind Countable Fluid bottle position.	COUNTABLE FLUID
Countable Fluid (bottle and tubing connector) : located on the cap of the Countable Fluid bottle and tubing connector on the instrument.	
Waste (lower): located on the work surface in front of the waste bottle position.	
Waste (upper): located on the vertical surface behind the waste bottle position.	
Waste (bottle and tubing connector): located on the cap of the waste bottle and tubing connector on the instrument.	

The Countable Computer

The computer that comes as part of your Countable System includes a desktop computer, monitor, keyboard, and mouse. The desktop computer comes ready to go with Windows 11 and the Countable Control Software already installed. Specifications for the Countable Computer are listed below:

Item	Description
Operating system	Windows 11
Processor	Intel i7
Memory	128 RAM
Storage	512 GB OS 2 x 4 TB SSD 22 TB HHD

A recommended layout for the computer is shown below:



The Countable Software

The Countable Control Software helps you operate the Countable System and set up experimental runs. The Countable Control Software also:

- Analyzes images and reports molecule counts.
- Measures the connectivity of targets across channels when the optional linkage analysis feature is enabled.
- Provides sample quality metrics and molecule counts in molecules per 50 μL.

Countable consumables, accessories, and reagents

You'll need Countable Fluid (KT0007) to operate the Countable System. Six (6) bottles of 500 mL Countable Fluid is included when you receive your Countable System. To purchase more, contact us at success@countablelabs.com.

The following Countable consumables and reagents will also be required to operate and analyze samples on your Countable System. Check out the Countable PCR Reaction Preparation User Guide for more information on using Countable consumables, accessories, and reagents.

	You should only use Countable consumables, accessories, and reagents to
Note	operate the Countable System. If other products are used, you'll need to
	validate their performance for your application.

Cat#	Name	Description
KT0001	Countable Matrix Consumables Kit	 Single-use consumables for matrix formation. Supports 192 reactions per kit and includes: Matrix Column Strips and Strip Seals Matrix Tube Strips and Strip Caps
KT0002	Countable Swing Bucket Kit	Swing bucket inserts compatible with Eppendorf 5430/5430R centrifuge fitted with S-24-11-AT swing bucket rotor. Each spin holds 48 samples (8 samples × 6 buckets).

Cat#	Name	Description
ктоооз	Countable Matrix Kit	Reagents for Countable Matrix formation. Includes MR01, MR02, MR03. Supports 192 reactions per kit.
KT0004	4X Countable PCR Mix Kit	Amplification mix for 50 µL reactions. Compatible with hydrolysis probes and Universal Multiplex chemistry. Supports 192 reactions per kit.
KT0005	Universal Multiplex Set A Kit	4-plex multiplexing chemistry using standard primers—no custom probes required. Designed for compatibility with 4X Countable PCR Mix. Supports 192 reactions per kit.
KT0007	Countable Fluid Kit	Maintenance fluid for routine use with the Countable System. One pack supports up to 1 month of operation.
ктооо9	Countable Control Assay Kit	Positive control assay kit for Countable PCR. Required for new user training or routine QC testing.

CHAPTER 5

Operating the Countable System

Turning on your Countable System

- 01 Turn on the Countable Computer by pushing the computer power button.
- 02 Log into the computer using your user credentials.
- 03 Power on your Countable Instrument by pressing the power button found on the right side of the instrument.



A LED indicator light at the bottom, right side of the front panel reflects the status of your instrument.

Countable System status	LED status
Power off	No lights
Power on/idle	Left light white
Running	All lights white
Error	Left light white, center and right light red

The table below summarizes different states for your Countable System.



04 Double-click on the Countable Control Software icon on your desktop to launch the application.



05 The Countable Control Software Main Menu will open.



Countable Control Software overview

The Countable Control Software Main Menu is used to operate and maintain your Countable System. It is a central hub to access the following core features of the system: Workbook maker, Image analysis, View results, Countable Fluid refill.

Tile	Function
Build Sample Workbook →	Build a Sample Workbook opens the Sample Workbook Maker that helps you create, organize, or import workbooks.
Image Samples	Image Samples starts imaging and analyzing your samples based on the experimental parameters defined in your workbook(s).
View Experiment → Reports	View Experiment Reports gives you access to the results of your analyses for your review. Processed data and outcomes are also stored here for easy access.
$\stackrel{\frown}{\underset{Fluid}{\overset{\bullet}}}$	Refill Countable Fluid manages liquid levels for the on-board Countable Fluid bottle and waste bottles.

Building a Sample Workbook

A Sample Workbook is created using the Countable Control Software and defines imaging and analysis parameters. A Sample Workbook must be created with each imaging and analysis run.

You can choose between two modes when building a Sample Workbook:

- **Position Mode**: The position of the Matrix Tube Strip in the carousel is fixed. You'll need to enter sample information based on the position of the Matrix Tube Strip in the carousel.
- Barcode Mode: The position of the Matrix Tube Strip in the carousel is flexible. You'll need to enter sample information and the corresponding barcode on the Matrix Tube Strip.

To build a Sample Workbook:

Ol Click on the Build Sample Workbook tile on the Countable Control Software Main Menu.



02 The Sample Workbook Maker Screen will open. If at any point you want to return to the Main Menu, click **Return to Main Menu**.

Countable Labs.		i – 🗆 ×
Sample Workbook Maker		E Import Countable Workbook
1. Select Mode Imaging Channels Imaging Channels Imaging Channel Details Imaging Channel Details Imaging Channels Imaging	Add Matrix Tube Strips	66
Return to Main Menu		SAVE

03 If you are creating a new Sample Workbook, skip to → <u>Defining Sample Workbook experimental</u> parameters

If you are importing a workbook that was created previously, select Import Sample Workbook.

04 Select the Sample Workbook you want to run in the dialog box that appears.

1. Select Mode	, Import Workbook	×	
Please select a mode	$\label{eq:constant} \leftarrow \ \rightarrow \ \lor \ \ \uparrow \ \ \ \ \ \ \ \ \ \ \ \ \ \ \$	Sample/Workbooks Y C Search Sample/Workbooks P	6 6
2. Select Imaging Channels	Organize • New folder	= · 🗆 🔮	
Ch01 Ch02 Ch03 Ch03 Ch03 Ch03 Ch03 Ch03 Ch03 Ch03	Countable Analysis Countable	Later monthes rype 244 5//(2025 835 AM Microsoft Everi Com.	
3. Enter Experiment Details Exp Name: Viser (Ontional):			
4. Use Training Samples (Optional)	> Whetwork		
5. Enable Linkage Analysis (Optional)	File name:	CSV File ("can) Concel Dpen Cancel	

05 Click Open.

Defining Sample Workbook experimental parameters

Countable Labs.		
Sample Workbook Maker		Import Countable Workboo
 Select Mode Select Imaging Choi Choi Select Mode Choi Select Mode Select Mode <li< th=""><th>Add Matrix Tube Strips</th><th>67 e7</th></li<>	Add Matrix Tube Strips	67 e7
Return to Main Menu		SAVE

_

OI Select **Position Mode** or **Barcode Mode** in the pull-down menu.

02 Select the imaging channels to use for imaging and analysis.

		import Countable workb
•		
1. Select Mode Position Mode *	Add Matrix Tube Strips	
	🕂 Add 2 - Strip(s) 🗙 Remove Last (🗙 Remove All	
2. Select Imaging Channels	Tube Position Sample Name	
✓ Ch01 ① ✓ Ch02 ① ✓ Ch03 ② ✓ Ch04 ①	Consumed Desiries of	^
	1 POI Sample()	Ť.
	2 POI_Sample02	
3. Enter Experiment Defails	3 PO1_SampleO3	
Exp Name: E2a	4 PO1_SampleO4	
User (Optional):	Carousel Position: 2	\$
	1 P02_SampleO5	
4. Use Training Samples (Optional) 🛛 🛈	2 P02_SampleO6	
	3 P02_SampleO7	
5. Enable Linkage Analysis (Optional) 🔵 0	4 P02_SampleO8	

Countable Labs.		i – 🗆
Sample Workbook Maker		Import Countable Workbo
1. Select Mode Position Mode • 2. Select Imaging Channels ⊘ Ch01 ⊘ Ch02 ♥ Ch03 ♥ Ch04	Add Matrix Tube Strips Add 1 - Strips) Remove Last Control Remove All Tube Position Sample Name	- -
3. Enter Experiment Details		
 Use Training Samples (Optional) Training Samples (Optional) Enable Linkage Analysis (Optional) 		
		_

03 Add a name for the experiment. You can also have the option to add user information.

04 Add the number of Matrix Tube Strips using the pull-down next to the 😌 icon.



- 05 Start entering sample details.
 - a In **Position Mode**, you must input sample names that correspond to the Matrix Tube Strip placed in a specific carousel position. You can also use default sample names.

Countable Labs.		i – 🗆 ×
Sample Workbook Maker		Import Countable Workbook
1. Select Mode Position Mode	Add Matrix Tube Strips	
2. Select Imaging Channels	Add 2 × Strip(s) Remove Last K Remove All Tube Position Sample Name	
	Carousel Position: 1	÷
3. Enter Experiment Details	2 P01_SampleO2 3 P01_SampleO3	
Exp Name: E2a User (Optional):	4 P01,Sample04 Carousel Position: 2	•
4. Use Training Samples (Optional)	1 P02_Sample05 2 P02_Sample06	
5. Enable Linkage Analysis (Optional) 🔘	3 P02_Sample07 4 P02_Sample08	
Return to Main Menu		SAVE

b In Barcode Mode you will input sample information and scan the corresponding Matrix Tube Strip barcode using a compatible USB barcode scanner. We recommend the Zebra Hand Held Barcode Reader (P/N LS2208).

mple Workbook Maker		Import Countable Workboo
1. Select Mode 🛛 Barcode Mode 🔹 😗	Add Matrix Tube Strips	
	🕂 Add 2 = Strip(s) 😣 Remove Last (😣 Remove All	
2. Select Imaging Channels	Tube Position Sample Name	
	Strip Barcode:	\$
 Enter Experiment Details 	2 P01_Sample02	
Exp Name: E2a	3 PO1_Sample03 4 PO1_Sample04	
User (Optional):	Strip Barcode:	:
	1 Power other a barcode	
4. Use Training Samples (Optional) 🛛 🔘 🕕	2 P02_SampleO6	
	3 P02_SampleO7 4 P02_SampleO8	
 Enable Linkage Analysis (Optional) 		

- 06 Optional: Select **Using Training Samples** if you are including a training sample to improve rare molecule detection. Refer to the Countable PCR Reaction Preparation User Guide for information on how to prepare this sample.
- 07 Optional: Select **Enable linkage analysis** to determine if two or more targets are physically connected in a single molecule to enable applications like measuring DNA integrity and fragmentation in a sample.
- 08 Click **Save** to create your Sample Workbook. You will be asked to give the workbook a name before returning to the Countable Control Software Main Menu.

> This PC > Local Disk (C) > CountableLabs > SampleWorkbooks	 ✓ C Search SampleWorkbooks 	×	
> This PC > Local Disk (C) > CountableLabs > SampleWorkbooks	✓ Ø Search SampleWorkbooks	0	
		~	
er		•	"o "o
ye Nome	Date modified Type		
	Working on it		
			· · · ·
20			÷ 0
ia) File (*.cov)			
	Save Cancel		
	2) Note	20 Exercised Statements (1994)	2) Anno Die monitor (1)00 Violarg un L.

Sample imaging

Starting sample imaging and analysis

To begin imaging and analyzing your sample:

01 Click on the Image samples tile on the Countable Control Software Main Menu.

Countable Labs.			i – 🗆 x
Let's star	t counting.		
Build Sample $ ightarrow$ Workbook	Image Samples $ ightarrow$	View Experiment \rightarrow Reports	$\stackrel{\text{Refill Countable}}{\text{Fluid}} \rightarrow$

02 The Imaging and Analysis Menu will appear. If at any point you want to return to the Main Menu, select **Return to Main Menu**.

Countable Labs.	i	- 🗆 ×
	Getting things started	
Countable Labs.	۵ Unice	k Door
Initialize	Make sure power is on and the door is closed. All hardware will be initialized, a checks will be performed.	nd system
Return To Main Menu	Select	Sample Workbook

O3 Click on the Initialize button to start the system initialization process. Confirmation will appear after the system has confirmed that the hardware is operating correctly, that there is enough Countable Fluid in the system, and that the Countable Computer has enough storage capacity to run your experiment.

Countable Labs.		i – 🗆 ×
HOME / ACQUISITION	Getting things started	
		🛱 Unlock Door
and the second second second	Initializing hardware	0
	Checking fluidics	0
Countable Labs.	Initializing sample handler	0
_	Checking disk integrity	0
Initialize	The system has been successfully initialized.	
Return to Main Menu		Select Sample Workbook

04 Click **Select Sample Workbook** to choose and load the workbook you want to use for the imaging and analysis.

HOME / ACQUISITION	Let's review t	he run	
Copen File	C. Open	×	
	← → ~ ↑ ► This PC → Local Disk (C3 → CountableLabs → SampleWorkbooks)	✓ C Search SampleWorkbooks ,P	
Sample Info	Organize * New folder	≡ · □ 0	rameters
	Videos / Name	Date modified Type Size	
	Desktop / Demo.csv	4/17/2025 5:04 PM Microsoft Excel Com	
	SampleWorkbooks		
	Mdeo 1		
	la loga		
	Logs		
	• This PC		
	> = user provider		
	The second		
	-		
	File name: Demo.cov	 CSV File (*.csv) 	
		Open Cancel	
		Optional I	nfo

05 Review the Run info, Imaging parameters, and Sample info parameters for your experiment in the Acquisition menu.

HOME / ACQL	IISITION			L	et's rev	view the	e run			
Open File	C:\Coun	tableLab	os\SampleV	Vorkbook	s\Demo.c	SV				
Sample Info									Imaging Paramet	ers
Carousel Position	Tube Position	Barcode	Sample Name	Channel01	Channel02	Channel03	Channel04		Channel Options	Channels to Scan
1	1		P01_SampleO1					^	Channel01	×
	2		P01_Sample02						ChannelO2	×
	3		P01_Sample03						Channel03	×
	4		P01_SampleO4						Channel04	×
2	1		PO2_SampleO5							
	2		P02_SampleO6						Run Info	
	3		P02_Sample07						Version	vl
	4		P02_Sample08						User	
3	1								Exp ID	E2a
	2								Mode (All/Slice)	All
	3									
	4								Optional Info	
4	1								Linkage Analysis ("x" or B	LANK)

- a If all the parameters look good, click **Confirm** to start sample loading.
- b If you want to edit some parameters, you'll need to go back to go back to the Workbook Maker, import the workbook, make the desired adjustments, and save the workbook.
- 06 The Loading Samples screen will appear. Click Load samples to unlock the Countable Instrument door.

	Let's load your samples
0.000	Sample Mapping
	Carousel Position Barcode Tube Positions Sample Names
	1 1,2,3,4 P01_Sample01,P01_Sample02,P01_Sample03,P01_Sample04
	2 1.2.3.4 PO2_SampleO5.PO2_SampleO6.PO2_SampleO7.PO2_SampleO8
Load Samples	

07 Open the Countable Instrument door and start loading your samples into the carousel, starting from position 01. If you have barcodes on your Matrix Tube Strips, make sure that they are facing left.



- 08 After all of your samples have been loaded, close the Countable Instrument door.
- O9 If you are operating in Barcode mode, the instrument will automatically start scanning the Matrix Tube Strip barcodes and will map the loaded strips. Confirmation that all barcodes were successfully scanned will appear.

Countable Labs.			I	-	×
	Let's load your sam	ples			
000	Sample Mapping				
	Carousel Position Barcode	Tube Positions	Sample Names		
	1 1JH13C	1,2,3,4	P01_SampleO1,P01_SampleO2,P01_SampleO3,P01_SampleO4		
A C A C	2 1JH133	1,2,3,4	P02_Sample05,P02_Sample06,P02_Sample07,P02_Sample08		
Load Samples					
Barcodes were successfully scanned and no issues were found.					
Return To Main Menu				Begin	

10 If there are barcode scan errors the strip with the issue will be highlighted in the Sample Mapping table. You'll need to check the mismatched Matrix Tube Strip(s) before continuing.

	Let's	s load your :	sampl	es				
0		Sample Mappir	g					
		Carousel Position	Barcode	Tube Positions	Sample Names			
	e la companya de la c	1	1JH13C	1,2,3,4	P01_Sample01,P01_Sam	ple02.P01_Sample03	3,P01_Sample	D4
	~	2		1,2,3,4	P02_Sample05,P02_Sam	nple06,P02_Sample	07,P02_Samp	\leftarrow
	8							
Load Samples								
There is a mismatch between the barcodes information in the sample sheet. Edit and r	scanned and eload the sample sheet							

11 Click **Begin** to start the imaging and analysis process.

Countable Labs.	i - 🗆	×
	Let's load your samples	
0 0	Sample Mapping	
AN AN AD	Carousel Position Barcode Tube Positions Sample Names	
19 0 2 0 2 0 0 M	1 1JH13C 1.2.3.4 P01_Sample01.P01_Sample02.P01_Sample03.P01_Sample04	
	2 1JH133 1.2.3.4 P02_Sample05.P02_Sample06.P02_Sample07.P02_Sample08	
Load Samples		
Barcodes were successfully scanned and no issues were four		
Return To Main Menu	Begin	\leftarrow

During sample imaging and analysis

The Countable Instrument door will remain locked while your samples are being imaged and analyzed.

01 The Acquisition Menu will show a progress bar with information about the status of your imaging and analysis, and provide a projection of when sample imaging will be completed.

HOME / ACQUISITIO	N		Makin	g it count						
Current Sample			P01_Sample04						û Unlock D)oor
Chantura			Dunning EO%	Sample Resu	ılts				~~~~~	🗹 Autosci
Status			Running 50%	Tube Position	Sample Name	Ch01 Counts	Ch02 Counts	Ch03 Counts	Ch04 Counts	
SAMPLES SCANNED TIME R	EMAINING (HRS:MINS)			Carousel Position	n:1 Barcos	de: 1JH13C				٥
4 /8 0:1	5		Abort	WellO1	P01_SampleO1	94020	80008	82629	116465	
				WellO2 WellO3	P01_SampleO2 P01_SampleO3	95173	107946	105963	93597	
Composite Images				WellO4	P01_Sample04		,	,	. J.	
ChUI	ChUZ	CHUS								

- 02 Snapshots of each sample will also appear for each of the previously selected channels.
- 03 The table below describes different fields that you'll see during the run.

Field	Description
Current sample	Indicates the sample that is currently being imaged.

Field	Description
Status	Indicates the status and percent completion of the analysis.
Samples scanned	Indicates the number of samples already scanned out of the total number of samples in the acquisition.
Time remaining	Estimates the remaining time of the imaging in (HRS:MIN).
Sample results	Displays the total counts per sample per 50 μL , for each channel selected.
Counts	Reports the total counts per 50 μL per channel for each sample.

04 If you need to stop the imaging and analysis mid-run, click **Abort**, then select **Return to Main Menu.** You'll be able to start a new imaging run from here.



At the end of the run

- 01 When the Countable System is done imaging and analyzing your samples, you'll be able to open the Countable Instrument door.
- 02 Select Unlock door and open the door.

HOME / ACQUI	SITION		Making it count								
Current Samp	le		P02_Sample08					-	Unlock D	Noor (
				Sample Resu	Sample Results				💼 💼 🛛 Autoscroli		
Status		Scanning C	ompleted 100%	Tube Position	Sample Name	Ch01 Counts	ChO2 Counts	Ch03 Counts	Ch04 Counts		
AMPLES SCANNED	TIME REMAINING (HRS-MINS)			Carcusel Positio	n:1 Barcoo	ie: 1JH13C				0	
8 /8	0:00		Abort	WellO1	P01_SampleO1	94020	80008	82629	116465		
				WellO2	P01_Sample02	117866	88640	85479	89867		
Composite Imag	jes			WellO3	P01_Sample03	95173	107946	105963	93597		
Ch01	Ch02	Ch03	ChO4	Well04	P01_SampleO4	109266	111921	80374	101991		
		0.100	onou -	Carousel Positio	n: 2 Barcoo	de: 1JH133				\$	
S. CARLEY	Conservation of the			WellO1	P02_SempleO5	86100	101028	111284	119412		
				WellO2	P02_SampleO6	84438	103965	114694	83667		
	1000			WellO3	P02_Sample07	109427	103497	100063	89131		
			2011 - 101 - 101 - 104 101 - 102 - 102 - 104	WellO4	PO2_SampleO8	93832	96010	113698	102641		
Service .	States H		201220								

- 03 Remove the Matrix Tube Strips from the Countable Instrument carousel.
- 04 Close the Countable Instrument door and click Return to Main Menu.

Viewing experiment reports

You can view a summary of the run and the data collected in the experiment report. The data can also be accessed as .html and .csv files. See \rightarrow <u>Output data files</u> for information on accessing the files in the [CountableLabs] folder.

To view a report in the Countable Control Software:



01 Click on the View experiment reports tile in the Countable Control Software Main Menu.

O2 Select the results you want to view in the selection screen.

n								□ ×
Cou	intable Labs.							Î
R	Count		reie Summo	*17				
oric	Summary	Quality Check Resul		' y				
	Experin	nent Overvie	W ?		QC Su	mmary 🕐		
	E	xperiment	E2a		No issues w	ere detected for this e	xperiment.	
		2025-04-10 12:	4:19.072271					
		Version	v0.1	v0.11				
	N	o. Samples	8					
	No	. Channels	4					
	Link	age Analysis	Off					
	Sample	Summary 🛛						
		Barcod	le Carousel Position	Tube Position	Sample Name	Training Sample	Flags	
		+ 1JH130	C 01	01	P01_SampleO1			
		+ 1JH130	C 01	02	P01_SampleO2			
		+ 1JH130	C 01	03	P01_Sample03			
		+ 1JH130	C 01	04	P01_Sample04			
								*

03 Once you've finished viewing your report, close the report and click Return to Main Menu.

Countable Analysis Summary

The Countable Analysis Summary can be viewed in the Countable Control Software or as an .html file. The summary contains the following sections:

- Summary
- Quality check
- Results

Summary

The Summary tab includes the following panels:

- **Experiment overview**: An overview of the analysis parameters, including the number of samples and channels imaged, and whether additional analysis like linkage analysis was turned on. Information in this section should match the Sample Workbook.
- **QC summary:** A list of any warnings that occurred during the run, and recommendations for how to address the warnings.
- Sample summary: A list of all the samples imaged. If Training Samples were included in the run, a colored square icon matching the channel selected will appear.

If a warning was flagged for a specific sample, the panel will default to an expanded view so you can get additional warning information. Warnings can be for an entire sample or can be channel specific.

The following flags will also appear in the panel so you will know how to approach your results:

Issues 🗢 : Reported counts for this sample/channel should not be used.

Observations **A** : Additional investigation of sample/channel quality should be performed before the reported counts are used.

LUDS								
Countab	ole Ar	nalysis Sum	nmary					
Summary Qua	ality Check	Results			00.5	Immory 2		
Lypenine					90.90			
Experiment			E2a		No issues	were detected for this ex	periment.	
	Date		04/15/2025	11:30				
User			v011					
No. Samples		12	12					
No. Channels		4						
L	.inkage Anal	lysis	On					
Sample S	Gumma	Carousel Position	Tube Position	Barcode	Sample Name	Training Sample	Flags	
Sample S	umma	Carousel Position	Tube Position	Barcode 1JH13C	Sample Name PO1_SampleO1	Training Sample	Flags	
Sample S	umma + +	Carousel Position 01 01	Tube Position 01 02	Barcode 1JH13C 1JH13C	Sample Name PO1_SampleO1 PO1_SampleO2	Training Sample	Flags	
Sample S	Summa	Carousel Position 01 01 01	Tube Position 01 02 03	Barcode 1JH13C 1JH13C 1JH13C	Sample Name PO1_SampleO1 PO1_SampleO2 PO1_SampleO3	Training Sample	Flags	
Sample S	6umma	Carousel Position 01 01 01 01 01	Tube Position 01 02 03 04	Barcode 1JH13C 1JH13C 1JH13C 1JH13C	Sample Name PO1_SampleO1 PO1_SampleO2 PO1_SampleO3 PO1_SampleO4	Training Sample	Flags	
Sample S	6umma + + + +	Carousel Position 01 01 01 01 01 01 01 02	Tube Position 01 02 03 04 01	Barcode 10H13C 10H13C 10H13C 10H13C 10H13C	Sample Name PO1_SampleO1 PO1_SampleO2 PO1_SampleO3 PO1_SampleO4 PO2_SampleO5	Training Sample	Flags	
Sample S	6umma + + + + + + +	Carousel Position 01 01 01 01 01 02 02	Tube Position 01 02 03 04 01 02	Barcode UH13C UH13C UH13C UH13C UH13C UH133 UH133	Sample Name POI_SampleO1 POI_SampleO2 POI_SampleO3 POI_SampleO4 PO2_SampleO5 PO2_SampleO6	Training Sample	Flags	
Sample S	Summa + + + + + + +	Carousel Position 01 01 01 01 02 02 02 02	Tube Position 01 02 03 04 01 01 02 03	Вагсоde	Sample Name POL_SampleO1 POL_SampleO2 POL_SampleO3 POL_SampleO4 PO2_SampleO5 PO2_SampleO6 PO2_SampleO6	Training Sample	Flags	

Quality Check

The Quality Check tab gives you fluorescent intensity of samples per channel in the Signal and Background panel. The panel displays histograms of the fluorescence intensity distribution of the target signal (colored) and background (grey), with the Intensity Distribution (ID) score appearing at the upper-left corner.

Countable	Labs.								
	Summary Quality Check		Summ	nary					
	Signal and Back	ground	?						
		Ch01 Score:	L.	ChO2	Score	Ch03	Score:	Ch04	
	PO1_SampleO1	94.9 3 4	5	nan 3 4	98.9	4	nan 5 3	4 5	
	P01_Sample02	99.3 3 4	5	99.7 3 4	Score nan 5 3	4	Score: nan	4 5	
	P01_Sample03	Score: 75.5		Score: 99.8	Score 99.4		Score: 99.6		
	P01_SampleO4	3 4 Score: nen	5	3 4 Score: 99.0	5 3 Score nan	4	5 3 Score: 99,9	4 5	
	P02_Sample05	3 4 Score: nan	5	3 4 Score: nan	5 3 Score 98.6	4	5 3 Score: 99.6	4 5	\bigcirc

Here are some guidelines on using the histograms to evaluate the sample and assay quality:

- 01 The signal and background distributions should be as distinct as possible.
- 02 Secondary peaks and shoulders between the signal and background intensities can indicate nonspecific primer/probe interactions. Normal intensity units for the signal should be greater than 4.
- 03 A long intensity tail connecting the signal and background distribution may indicate spatial variations, such as PCR inhibition.

Results

The Results tab uses a bar chart to visualize the counts per 50 μ L for each sample. You can toggle each channel between linear- and log-scale views depending on your preference.

Countable Labs.							
	Countable An	alysis Sumr	mary				
	Counts per 50 µ	L ?					
		Ch01	Ch02 Linear	Ch03	Ch04		
	P01_Sample01	560	0	9662	1		
	P01_Sample02 93	538	9346	1	0		
	P01_Sample03 1329		9590	9635	9568		
	P01_SempleO4 15		9551	n	9538		
	P02_Sample05 1		1	9985	9935		
	P02_SampleO6	781	4	5	9815		
	P02_Sample07 153		9577	9633	0		
	P02_SampleO8	947	9942	9927	0		
	P03_SampleO9 821		809	846	869		
	P03_Sample10 93	544	9639	9674	9664		
	P03_Sample11	95998	96386	95987	97362		
	P03_Sample12 15		36	9	79		

Output data files

Several data files are generated when your analysis is done. Data files can be found in the [CountableLabs] folder, which contains subfolders organized by the Countable Control Software versions. Each version folder contains data and metadata folder which each contains folders that correspond to the unique experiment IDs you designated.

Folder	Description
Data	Contains results including the Countable Analysis Summary (.html), the Countable Data Report (.csv) and Snapshot folder. It also contains the Linkage Summary (.csv) if the linkage option was selected for the experiment.
Metadata	Contains folders (.zip) for each sample of the experiment. Each sample contains files generated during the analysis process. It also contains one folder (.zip) for experiment. These folders can be shared with Countable Labs if support and troubleshooting is needed.

Folder and file names for each experiment includes the Experiment ID, date, and time using the following format: ID_YYMMDD_HHmm. For example. "Expla_250409_1125" is experiment Expla that was run on April 9, 2025 at 11:25 am.

Countable Data Report

The Countable Data Report (.csv) contains the sample counts, in counts per 50 μ L, in a format that can be opened in a spreadsheet software. The report can be accessed from the "View experiment results" tile of the Countable Control Software and in the [CountableLabs] folder.

Countable Data Report file names follow this format: ExperimentID_timeStamp_ CountableDataSummary.csv. Each line corresponds to one sample. An example of a Data Report is captured below.

CarouselPos	TubePosition	StripBarcode	SampleName	SampleID	Ch01_Count	Ch02_Count	Ch03_Count	Ch04_Count	Ch01_IDsco	or Ch02	_IDscor	Ch03_ID	scor C	h04_IDsco	Ch01_Signal	Ch02_Signal	Ch03_Signal	Ch04_Signal	Ch01_Backg	Ch02_Backgi	Ch03_Backg	Ch04_Backg
1	1	1JH13C	P01_Sample	1JH13C_01	9619	0	9646	0	94.1	1 NA			99 N	A	13213	0	20893	0	563	39	609	201
1	2	1JH13C	P01_Sample	1JH13C_02	9346	9355	1	0	99.3	3	99.7	NA	N	A	16376	21057	5764	0	352	336	3	3
1	3	1JH13C	P01_Sample	1JH13C_03	1317	9577	9622	9555	NA		99.8	ę	9.4	99.6	1653	21597	25509	52179	82	630	619	277
1	4	1JH13C	P01_Sample	1JH13C_04	45	9547	5	9533	NA		99	NA		99.8	1757	24470	22902	57201	66	590	4	300
2	1	1JH133	P02_Sample	1JH133_01	0	1	9982	9942	NA	NA		ę	98.6	99.7	0	14570	24752	50309	47	31	603	203
2	2	1JH133	P02_Sample	1JH133_02	9765	5	3	9808	91.6	6 NA		NA		99	20024	44762	17629	59716	624	32	4	319
2	3	1JH133	P02_Sample	1JH133_03	614	9563	9618	0	NA		98.4	ş	98.2 N	A	1735	20548	23342	0	75	569	592	15
		4111400	D00 0 1 1	4111400 04	0054	00.45	0007								00070	05443	00700			100		

The table below summarizes the columns included in the Data Report.

Field	Description
CarouselPosition	Position of the Matrix Tube Strip on the carousel.
TubePosition	Position of the sample within the Matrix Tube Strip.
StripBarcode	Strip Barcodes on Matrix Tube Strips.
SampleName	Name of the sample provided in the Sample Workbook.
SampleID	Unique sample identifier (StripBarcode_TubePosition).
CountsPer50µL	Number of targets per 50 μ L detected in each channel.

Field	Description
IDscore	Intensity Distribution (ID) score: A unitless metric with a scale from $0-100$ that describes the distribution of fluorescence intensity of the detected single molecules per channel. Reported for samples with >3000 and <1 M counts. The higher the metric, the better the signal distribution is from the background. The metric is influenced by a combination of assay and sample quality. The optimal score should be >90.
Signal	Median of fluorescence signal intensity from single molecules per channel.
Background	Average fluorescence intensity of the rest of the tube with no identified single molecule signals per channel.
lssues	Issue with a sample preventing the generation of reliable counts per channel.
Observations	Observed warning signs in a sample which may lower accuracy of the counts in specified channels.

Countable Snapshots

Countable Snapshots are images of the cross section at the center of the sample tube per selected channel, taken by the Countable Instrument during imaging and analysis. They can be used to visually inspect the sample tubes for troubleshooting purposes.

The Experiment ID, Sample ID, and Sample Name can be found on the top left corner of each snapshot image. You can also find the name of each channel at the bottom of the image. Countable Snapshot files can be found in the [CountableLabs] folder and can be accessed from the "View experiment reports" tile of the Countable Control Software Main Menu.



Linkage Summary

The Linkage Summary (.csv) is created when you enable the linkage analysis option when building a Sample Workbook. The file contains molecule counts per 50 μ L for all possible channel configurations (i.e., all individual channels and their combinations) for each imaged sample. Linkage Summary files can be found in the [CountableLabs] folder and can be accessed from the "View experiment reports" tile of the Countable Control Software Main Menu.

Linkage Summary file names use the following format: ExperimentID_timeStamp_ CountableLinkageSummary.csv.

The Linkage Summary contains the following information:

- **Sample information**: Carousel position, Matrix Tube Strip position, Sample Name, and Sample ID.
- Molecular counts per 50 μL for all channel configurations: For each sample, the listed counts relate to the corresponding configuration. For example:
 - Ch01_only: counts from Ch01 only (no linkage).
 - Ch01_02: counts from Ch01 and Ch02, only in this configuration.
 - The sum of counts in (Ch01_only, Ch01_02, Ch01_03, Ch01_04, Ch01_02_03, Ch01_02_04, Ch01_03_04, and Ch01_02_03_04): the total counts in Ch01.
 - The sum of counts from all configurations: the total count per sample.
- **Flags**: Warning that there was a failure in the linkage analysis process.

CarouselPos TubePositio	on StripBarcod	SampleNam	e SampleID	Ch01_only	Ch02_only	Ch03_only	Ch04_only	Ch01_02	Ch01_03	Ch01_04	Ch02_03	Ch02_04	Ch03_04	Ch01_02_03	Ch01_02_04	Ch01_03_04	Ch02_03_04	Ch01_02_03	Flags
1	2 1JH13C	P01_Sample	(1JH13C_02	28	37	0	0	9317	C	0	(0 0	0	1	0	0	0 0	0	
1	3 1JH13C	P01_Sample	(1JH13C_03	30	72	89	37	8	3	8	74	1 57	87	7	4	7	8105	1250	
1	4 1JH13C	P01_Sample	(1JH13C_04	9	96	2	84	2	C	C	:	l 9413	1	0	34	0) 1	0	
2	1 1JH133	P02_Sample	(1JH133_01	0	0	82	42	0	0	C	(0 0	9899	0	0	0) 1	0	
2	2 1JH133	P02_Sample	(1JH133_02	161	0	0	203	2	C	9602	. (0 0	0	0	0	0) 3	0	
2	3 1JH133	P02_Sample	(1JH133_03	7	115	146	0	3	27	' C	8868	3 0	0	577	0	0	0 0	0	
2	4 1JH133	P02_Sample	(1JH133_04	133	97	159	0	216	136	i 0	166	6 0	0	9466	0	0	0 0	0	
3	1 1JH134	P03_Sample	(1JH134_01	7	0	6	21	0	C	0	() 5	3	3	5	38	3 28	762	
3	2 1JH134	P03_Sample	1JH134_02	102	32	63	76	42	20	9	12	2 19	24	47	34	52	411	9041	
2	3 110134	P03 Sample	110124 02	405	477	400	1422	242	170	171	100	277	206	424	500	271	706	02072	

Metadata

During data processing, a metadata folder is automatically generated. This folder contains more information about the data processing for a specific scan and analysis. Share the contents of this folder with Countable Labs support when you need troubleshooting help.

Refilling Countable Fluid

Matrix Tube Strips are immersed in Countable Fluid when they are imaged in the imaging chamber by the Countable Instrument. The Countable Control Software automatically tracks how much Countable Fluid you need for each experiment and will let you know if you need to replace the Countable Fluid bottle.

Caution

The Countable System damage

The waste bottle must be emptied every time the Countable Fluid is refilled.

To Refill Countable Fluid:

01 Click on the Refill Countable Fluid tile in the Countable Control Software Main Menu.



- i 🗆 × Countable Labs. Countable Fluid Bottle Exchange Imaging Fluid Bottle
 O Waste Fluid Bottle û Unlock Door ✓ 1. Unlock the door and manually lift the instrument cover to open. 2. Push tab and lift connector to unclip the empty Countable Fluid bottle. 3. Remove the empty Countable Fluid bottle from the instrument. ✓ 4. Swap the cap from the empty bottle to the full bottle. 5. Place the full Countable Fluid bottle in the Countable Fluid station. 6. To reconnect, press Countable Fluid bottle connector until it clicks. 7. Manually close the instrument cover. Manually lift the instrument cover to open Test Fluidics • • • Return To Main Menu
- 02 A screen will appear with a video that will take you through the steps listed below on how to refill the Countable Fluid.

- 03 Open the Countable Instrument door to reveal the amber Countable Fluid bottle and white/clear waste bottle.
- 04 Press the silver notch on the Countable Fluid bottle to engage it. You will feel or hear a click.
- 05 Gently remove the male connector. There should be minimal pressure if the notch is correctly engaged.
- 06 Remove the empty bottle from the Countable Instrument.
- 07 Unscrew and remove the cap on a new bottle of the Countable Fluid (KT0007).



- 08 Unscrew and remove the cap from the empty Countable Fluid bottle. Watch out for any residual fluid on the internal tubing.
- 09 Transfer the cap from the empty Countable Fluid bottle to the new Countable Fluid bottle.
- 10 Place the new Countable Fluid bottle back onto the Countable Instrument in the correct bottle position.
- 11 While supporting the Countable Fluid bottle, connect the tubing connector on the instrument to the connector on the bottle cap. Listen or feel for a click to ensure the connectors are fully closed to prevent leaks.
- 12 Close the Countable Instrument door.
- 13 Click Return To Main Menu.
- 14 Cap the empty Countable Fluid bottle and properly dispose of it according to your lab protocols in the proper waste container.

Emptying the waste bottle

The waste bottle contains used Countable Fluid and is located next to the Countable Fluid bottle on the Countable Instrument.

Caution



The waste bottle must be emptied every time the Countable Fluid is refilled.

To empty the waste bottle:

- Ol Click on the Refill Countable Fluid tile in the Countable Control Software Main Menu.
- 02 Open the Countable Instrument door to reveal the amber Countable Fluid bottle and white/clear waste bottle.
- 03 Press the silver notch on the waste bottle to engage it. You will feel or hear a click.
- 04 Gently remove the male connector. There should be minimal pressure if the notch is correctly engaged.
- 05 Remove the waste bottle from the instrument.
- 06 Unscrew the cap and dispose of the used Countable Fluid. The Countable Fluid waste can be poured into a larger fluid waste container for chemical disposal according to your lab protocol and local guidelines.
- 07 Screw the cap onto the waste bottle.
- 08 Place the empty waste bottle with cap back onto the instrument in the correct bottle position.
- 09 While supporting the waste bottle, connect the tubing connector on the instrument to the connector on the bottle cap. Listen or feel for a click to ensure the connectors are fully closed to prevent leaks.
- 10 Close the Countable Instrument door.
- 11 Click Return To Main Menu.

CHAPTER 6

Maintenance

Cleaning the Countable System

Tu

Risk of personal injury and the Countable System damage

Turn off the Countable Instrument and the Countable Computer before cleaning the instrument.

Once a month:

Caution

01 Apply 70% isopropyl alcohol to a lint-free wipe.



02 Wipe down the Countable covers and working surfaces.



Cleaning carousel spills

Samples might contain Countable Fluid remnants when they are transferred from the imaging chamber to the carousel. The remnants may accumulate over time.

To clean the carousel:

- 01 Remove the carousel.
- 02 Clean it with 70% isopropyl alcohol on a lint-free wipe.
- 03 Replace the carousel in the instrument.



Risk of personal injury and the Countable System damage

- Apply cleaning solution to the wipe and not directly to the instrument to avoid excess solution leaking into the instrument and damaging components.
- Do NOT use bleach to clean the Countable System.

Removing equipment from use for repair or disposal

Only Countable Labs personnel should prepare the system for removal from use, transportation, or disposal. Reach out to Countable Labs for information on how to prepare for their visit.

CHAPTER 7

Troubleshooting

Every now and then, you may need to troubleshoot your Countable System. Some common fixes are listed in this chapter. If you still need help, contact Countable Labs support at success@countablelabs.com.

Sending data to Countable Labs for troubleshooting or support

If you're asked to send data files from your Countable Computer as part of the support process, send requested files to success@countablelabs.com. If requested, information about your instrument serial number can be found on a label at the back of the instrument.

System initialization errors

Error	User action
Software crashes or doesn't initialize	 Restart the computer and try to initialize the system again. If the issue persists, contact Countable Labs Technical support at <u>success@countablelabs.com</u>.
Unable to connect to the camera	01 Unplug and replug the blue USB connection from the desktop computer to check the camera USB connection.02 Restart the Countable Control Software.

Error	User action
Unable to connect to barcode scanner	01 Unplug and replug the USB from the desktop computer to check the barcode scanner USB connection.02 Restart the Countable Controls Software.
One or more lasers are not warmed up	O1 Wait until the LED indicator lights are no longer flashing (~1 minute).O2 Try to initialize the system again.

Workbook is not accepted

Description	User action
A workbook is not accepted by the Countable Control Software	O1 Make sure that the Matrix Tube Strips are placed in the appropriate position and orientation in the carousel.O2 Make sure that the samples listed in the workbook match the samples placed in the carousel.

Fluidics

The Countable System cannot start a new acquisition if the Countable Fluid bottle or waste bottle is not properly connected. The Countable Control Software will send an alert in the Countable Control Software if there is an issue.

Description	User action							
The tubing connector have not properly engaged the caps of the Countable Fluid bottle and/or waste bottle	01 Remove and reattach connectors again.02 You should feel or hear a small click when the connectors are engaged.03 Rerun the instrument initialization.							
The Countable Fluid bottle may not have enough fluid to fill the internal Imaging Chamber	 Check if the Countable Refill bottle is near empty. Refill the Countable Fluid bottle, following steps in Chapter 5. Empty the waste bottle. Rerun the instrument initialization. 							
Persistent error messages	 Contact Countable Labs Technical Support at success@countablelabs.com if: 01 The system thinks the Countable Fluid or waste bottle is missing when they are both present. 02 The system thinks the Countable Fluid bottle is empty when they are full. 03 The system incorrectly flags a full waste bottle. 							

Countable Analysis Summary issues and observations

The tables below summarize issues and observations you may see in the Countable Analysis Summary (.html) file.

Reported issue	Description	User action			
Analysis failed	The Countable System analysis did not run to completion.	Contact Countable Labs support to troubleshoot.			
Low Intensity Distribution (ID) Score	Intensity Distribution Score is low.	Assay may require further optimization.			
Reported observation	Description	User action			
ChX high count accuracy alert (>1 M)	Detected counts are greater than 1 million and counting accuracy may be suboptimal.	Repeat experiment with less input DNA template to lower count.			
The Countable Control Software doesn't detect the Training Sample or asks you to provide an appropriate Training Sample for ChX to increase counting accuracy	The count of molecules per 50 μ L are less than 1000 counts and no training sample was provided OR the provided training sample had less than 1000 counts.	Repeat the experiment with a selected Training Sample included in the Sample Workbook. Refer to the Countable PCR Reaction Preparation User Guide for more details.			
Low signal	Fluorescence intensity of target DNA is low.	Assay may require further optimization.			
High background	Fluorescence intensity of matrix background is very high.	Assay may require further optimization.			



Countable Labs, Inc. 1810 Embarcadero Rd, Suite 200 Palo Alto, CA 94303

countablelabs.com