

## Countable Fluid Kit (KT0007)

For Research Use Only. Not for use in diagnostic procedures.

---

### Description

The Countable Fluid Kit (KT0007) includes six bottles of Countable Fluid used to operate the Countable System. Countable Fluid is required to fill the internal Imaging Chamber during sample imaging, enabling optical clarity and consistent refractive index for high-resolution 3D light sheet microscopy.

Each bottle contains pre-filtered imaging-grade fluid and is designed for seamless integration with the Countable System's fluidics module. Countable Fluid is automatically tracked by the Countable Control Software, which notifies users when refilling is required.

### Kit Contents

Component	Catalog #	Quantity
Countable Fluid	PR0014	6 × 500 mL

Supports up to a month of operations.

### Storage and Stability

- Store at ambient temperature (18°C to 25°C)
- Keep away from direct sunlight and other heat sources
- Do not freeze
- Shelf life: 2 years from date of manufacture

### Handling Instructions

None

### Workflow Overview

See Countable System Instructions for Use

## Safety Information

This kit contains flammable hazardous materials. Refer to the Safety Data Sheet. Avoid contact with skin and eyes. Avoid inhalation of vapor or mist. Wash hands with soap and water after handling. Handle at room temperature and use in accordance with standard lab safety practices.

## Additional Information

For a usage and setup, refer to the Countable System Instructions for Use.

For technical support, contact [success@countablelabs.com](mailto:success@countablelabs.com).

For order assistance, contact [order@countablelabs.com](mailto:order@countablelabs.com).

Visit [www.countablelabs.com](http://www.countablelabs.com)